RESCUE SAM

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Defibrillator

RESCUE SAM is a semi-automatic external defibrillator designed for providing simple and direct voice prompts and indications for a straightforward rescue operation. It is lightweight and battery powered for maximum portability.

RESCUE SAM is designed to treat Ventricular Fibrillation (VF) and Fast Ventricular Tachycardia (VT).

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SHOCK ENERGY ACCURACY:

According to IEC EN60601-2-4

SHOCK-TO-SHOCK CYCLE TIME: typically less than 20 seconds

PROTOCOL: voice prompts and indicators guide user through protocol

CHARGE CONTROL: automatic by software (arrhythmia detection
system and charging control)

CHARGE TIME FROM "SHOCK ADVISED":

<8 seconds, typical

VOICE INSTRUCTIONS: detailed voice messages guide responder
through use of the defibrillator

CONTROLS: shock button, info-button, On/Off button

INDICATORS: LEDs (different colors), info-button

CHARGE COMPLETE INDICATOR:

- Text prompt (press the shock red button)

- Flashing red shock button

DISARM: once charged, the RESCUE SAM disarms if:

- Patient's heart rhythm changes to non-shockable rhythm, or

- The SHOCK button is not pressed within 15 seconds after the is armed,
or

- The Power ON/OFF button is pressed to turn OFF, or

- The defibrillator pads are removed from the patient or the pads
connector is disconnected

SHOCK DELIVERY: shock is delivered if the SHOCK button is pressed
while the defibrillator is armed

SHOCK DELIVERY VECTOR: via adult defibrillator pads in the anterioranterior
(Lead II) position or via reduced-energy pediatric pads in the
anterior-posterior position

PATIENT ISOLATION: type BF

BATTERY: 12 Volt DC, 4.2 Ah, lithium manganese dioxide, disposable
long-life

BATTERY CAPACITY: min.200 shocks or 4 hours of operating time
continuous

AUTOMATED AND USER-ACTIVATED SELF-TESTS:

- Daily Self-tests: tests internal circuitry, waveform delivery system,
battery capacity, and software.

- Battery Insertion Test: upon battery insertion, extensive automatic selftests
and user-interactive test check device readiness.

DATA RECORDING AND TRANSMISSION:

- Optional special battery pack for data transmission

- Data stored: first 60 minutes of ECG and the entire incident's events
and analysis decision.

SIZE: 80x300x270 mm

WEIGHT: 2.2 kg with battery pack installed.

Standard set:

- main unit

- disposable defi pads adult

- disposable battery pack

Optionals:

- preconnected defi pads adult/pediatric

- disposable defi pads pediatric

- carrying bag

- wall mounted cabinet

- software with USB cable adapter for PC data
management

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